

ST1 Stretcher



Basic patient transport.

Stryker's ST1 is your basic transport stretcher. Designed to satisfy your fundamental department needs, this stretcher includes fifth wheel mobility for easier steering and Stryker's patented glideaway siderails.

Standard features

- Fifth wheel mobility
- Glideaway siderails
- 4 integrated IV receptacles
- Central lock mechanism

Optional features

- Basic mattresses: 8 cm x 62 cm x 193 cm (3" x 24" x 75.5") or 10 cm x 62 cm x 193 cm (4" x 24" x 75.5")
- Pop-up push handle options
- Restraint straps
- Hood base sticker (red, navy blue, gray)
- Upright oxygen bottle holder
- Defibrillator tray with chartholder
- Paper roll holder

Specifications

Model number

ST1	6300-000-100
-----	--------------

Overall length

217 cm (85.4")

Minimum height

56 cm (22")

Overall width

78 cm (30.7")

Siderails up	78 cm (30.7")
--------------	---------------

Siderails down	73.5 cm (28.9")
----------------	-----------------

Safe working load (SWL)

250 kg (551 lbs.)

Trend/Reverse Trend.

±16°

Patient surface

62 ± 0.5 cm x 193 ± 0.5 cm (24.4" x 76")
--

Warranty

Warranty may vary by country. Contact your local Stryker Medical representative for additional information.

Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, ST1 and ST1-X. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



PHMISC-Com-3_21867EN Rev 2
SDL 08/2020
2019-21867

Copyright © 2020 Stryker
stryker.com